**Alphatec Solus™ Anterior Lumbar Interbody Fusion (ALIF)
System**

GENERAL INFORMATION:

The Alphatec Solus Anterior Lumbar Interbody Fusion (ALIF) System is an intervertebral body fixation system consisting of implants with various heights and lordosis to accommodate individual patient pathology. System implants are manufactured from implant grade polyetheretherketone (PEEK Optima LT1), titanium (Ti-6Al-4V ELI) anchoring blades and tantalum radiographic markers. System instruments are manufactured from stainless steel. The Alphatec Solus implant is intended for use with supplemental spinal fixation. Specifically, the Alphatec Solus implant is to be used with the Alphatec's Zodiac® Spinal Fixation System, Aspida™ Anterior Lumbar Plating System, ILLICO® MIS Posterior Fixation System, ILLICO® FS Facet Fixation System, or the BridgePoint™ Spinous Process Fixation System.

WARNINGS:

1. The implants are provided sterile
 - a. Do not re-sterilize implants
 - b. Do not use implants if package is opened or damaged or if expiration date has passed.
2. Do not reuse implants or distal deployment tip.
3. Instruments must be cleaned and sterilized before each use per the parameters provided in the *STERILIZATION / RESTERILIZATION* section of this IFU.
4. Implant blades must be fully deployed into both vertebrae when implanted.
5. Potential risks identified with the use of this intervertebral body fusion device, which may require additional surgery, include: device component fracture, loss of fixation, pseudarthrosis (i.e., non-union), fracture of the vertebra, neurological injury, and vascular or visceral injury.
6. Implantation and positional adjustment of implants must only be performed using the special equipment and instruments specific to these devices. They must not be used with other instrumentation unless specifically recommended by Alphatec Spine Inc.
7. Do not use with excessively dense vertebrae. In patients where only one of the two vertebrae is excessively dense the safety feature built into the deployment wrench tip will not be effective.

PRECAUTIONS

1. Based on the dynamic testing results, the physician should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc., which may impact on the performance of the intervertebral body fusion device.
2. The implantation of the intervertebral body fusion device should be performed only by experienced spinal surgeons with specific training in the use of this device because this is a technically demanding procedure presenting a risk of serious injury to the patient.
3. The selection of the proper size, height and lordosis of the implant for each patient is crucial to the success of the procedure.
4. The Alphatec Solus Anterior Lumbar Interbody Fusion (ALIF) System has not been evaluated for safety and compatibility in the Magnetic Resonance environment.
5. Use fluoroscopic imaging when placing supplemental spinal fixation to verify there is no interference with the Alphatec Solus device.
6. Autograft must be used in conjunction with the Alphatec Solus Anterior Lumbar Interbody Fusion (ALIF) Spinal Spacer System to augment stability. Autograft should be packed inside the device prior to insertion, and around the device after insertion. The graft should extend from the upper vertebra being fused to the lower vertebra being fused. Use of this product without a bone graft may not be successful.
7. Bone density, quality, and/or stability may have adverse affects on the anchoring blades ability to be deployed or maintain purchase with the bone.
8. The Alphatec Solus Anterior Lumbar Interbody Fusion implants have not been evaluated for safety and compatibility in the MR environment. The Alphatec Solus Anterior Lumbar Interbody Fusion implants have not been tested for heating or migration in the MR environment.

INDICATIONS:

The Alphatec Solus Anterior Lumbar Interbody Fusion (ALIF) System is indicated for spinal fusion procedures in skeletally mature patients with degenerative disc disease (DDD) at one or two contiguous levels from L2 to S1. These DDD Patients may also have up to Grade 1 spondylolisthesis or retrolisthesis at the involved level(s). DDD is defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies. These patients should be skeletally mature and have had six months of non-operative treatment. The Alphatec Solus implant is intended to be used with autograft. The device is intended for use with supplemental fixation that is in addition to the integrated blades.

CONTRAINDICATIONS:

The Alphatec Solus Anterior Lumbar Interbody Fusion (ALIF) Spinal Spacer System is contraindicated for:

1. Use in the cervical spine.
2. Patients with osteopenia, bone absorption, bone and/or joint disease, deficient soft tissue at the wound site or probable metal and/or coating intolerance.
3. Patients with infection, inflammation, fever, tumors, elevated white blood count, obesity, pregnancy, mental illness and other medical conditions, which would prohibit beneficial surgical outcome.
4. Spinal surgery cases that do not require bone grafting and/or spinal fusion.
5. Patients resistant to following post-operative restrictions on movement especially in athletic and occupational activities.
6. Implantation using instruments from other systems.
7. Reuse, or multiple use of single use implant or distal tip instrument.
8. Prior fusion at the level to be treated.
9. Any condition not described in the indications for use.

POSSIBLE ADVERSE EFFECTS:

The following complications and adverse reactions have been shown to occur with the use of similar spinal instrumentation. These effects and any other known by the surgeon must be discussed with the patient preoperatively.

1. Initial or delayed loosening, bending, dislocation and/or breakage of device components.
2. Physiological reaction to implant devices due to foreign body intolerance including inflammation, local tissue reaction, and possible tumor formation.
3. Loss of desired spinal curvature, spinal correction and/or a gain or loss in height.
4. Infection and/or hemorrhaging.
5. Bone graft, vertebral body and/or sacral fracture, and/or discontinued growth of fused bone at, above and/or below the surgery level.
6. Non-union and/or pseudarthrosis.
7. Neurological disorder, pain and/or abnormal sensations caused by improper placement of the device, and/or instruments.
8. Scar tissue formation possibly causing neurological and/or vascular compromise.
9. Bone loss and/or decrease in density due to stress shielding.
10. Subsidence of the device into the vertebral body.
11. Revision surgery.
12. Death.

PREOPERATIVE MANAGEMENT

1. The surgeon should have a complete understanding of the surgical technique and of the device, indications, contraindications and applications.
2. The surgeon should have a complete understanding of the function and limitations of each implant and instrument.
3. All device components should be received and accepted only in packages that have not been damaged or tampered with. Components must be carefully handled and stored in a manner that prevents scratches, damage, and corrosion.
4. The condition of all implants and instruments should be checked prior to use. Damaged and/or worn implants or instruments should not be used.

INTRAOPERATIVE MANAGEMENT:

1. The surgical technique manual should be followed carefully.

- To prevent possible nerve damage and associated disorders, extreme caution should be taken to avoid the spinal cord and nerve roots at all times, especially upon insertion.
- Careful use of the implants and instruments should be taken. Misuse of the components could cause injury to the patient or operating personnel.

POSTOPERATIVE MANAGEMENT:

- The patient should have a complete understanding of and compliance with the purpose and limitations of the implant device. The surgeon should instruct the patient on how to compensate for any loss in range of spinal motion due to bone fusion.
- The surgeon should instruct the patient regarding amount and time frame after surgery of any weight bearing activity. The increased risk of bending, dislocation, and /or breakage of the implant devices, as well as an undesired surgical result are consequences of any type of early or excessive weight bearing, vibration motion, fall, jolts or other movements preventing proper healing and/or fusion development.
- In the case of delayed, mal-union, or non-union of bone, the patient must continue to be immobilized in order to prevent bending, dislocation, or breakage of the implant device. Immobilization should continue until a complete bone fusion mass has been developed and confirmed.
- Postoperative patients should be instructed not to smoke, consume alcohol, or consume non-steroidals and aspirin, as determined by the surgeon. Complete postoperative management to maintain the desired result should also follow implant surgery.
- Implanted devices should be revised or removed immediately if appropriate, upon a case of a non-union, pseudarthrosis or if the devices are bent, dislocated or broken.
- Retrieved implants should be properly disposed of in accordance with hospital procedures.

INSTRUMENT PREPARATION:

- Cleaning, inspection, lubrication, and sterilization must be performed by hospital personnel trained in the general procedures involving contaminant removal.
- Instruments must be cleaned prior to lubrication and sterilization.
- All instrument hinged, rotating, and articulating parts must be lubricated prior to sterilization with a water soluble and sterilizable lubricant intended for surgical instruments (Hinge-Free® for example).

CLEANING:

All instruments must be free of packaging material and bio-contaminants prior to sterilization. Cleaning, maintenance and mechanical inspection must be performed by hospital personnel trained in the general procedures involving contaminant removal.

For complex devices, such as those with, cannulas, hinges, retractable features, mated surfaces, and textured surface finishes, require special attention during cleaning. Manual pre-cleaning of such device features is required before automated cleaning processing.

Assure devices are in the fully extended, open position throughout cleaning. Certain instruments may require dismantling before cleaning.

All solutions for cleaning should be prepared per the manufacturer's recommendations.

Handle all products with care. Mishandling may lead to damage and possible improper functioning of the device.

Manual Cleaning Steps for Instruments

Step 1	Rinse devices in ambient temperature tap water to remove excess soil.
Step 2	Submerge device in enzyme solution, such as Polystica® 2X Enzymatic, and soak for a minimum of 10 minutes.
Step 3	Actuate and scrub the device using a soft bristled brush for a minimum of 2 minutes. If needed, actuate at several locations to access all surfaces. Use of a syringe (minimum of 50 ml) or water jet is recommended for hard to reach areas and repeat 3 times.
Step 4	Remove the devices from the detergent and rinse the

	devices in purified (RO/DI) water for a minimum of 1 minute.
Step 5	Submerge and actuate devices in cleaning solution such as Polystica® 2X Alkaline and sonicate for a minimum of 10 minutes
Step 6	Thoroughly rinse devices with RO/DI water to remove all detergent residues.
Step 7	Dry devices with clean, lint free cloth or filtered compressed air.

Automatic Washer Cleaning Steps

Step 1	Follow steps 1 through 5 of the <i>Manual Cleaning Steps for Instruments</i>
Step 2	Thoroughly rinse devices in ambient temperature tap water to remove detergent residuals.
Step 3	Place devices in fully extended open position into washer and process through a standard washer/disinfector instrument cycle.
Step 4	PreWash, cold tap water, 2 minutes.
Step 5	Enzyme wash (such as Polystica® 2X enzymatic), hot tap water, 1 minute
Step 6	Detergent wash (such as Polystica® 2X Alkaline), Hot tap water (66°C/150°F), 2 minutes
Step 7	Rinse 2x, hot tap water, 15 seconds
Step 8	Purified Water rinse, Hot tap water (66°C/150°F), 10 seconds
Step 9	Hot Air Dry, (115°C/239°F), for a minimum of 10 minutes

STERILIZATION / RESTERILIZATION:

All instruments are provided non-sterile and must be cleaned and sterilized before use. Instruments should be autoclave sterilized using the cycle parameters in the table below.

Alphatec instruments have been validated to achieve sterility using FDA cleared sterilization accessories (container and filters). FDA cleared filters should be used to achieve and maintain sterility after processing.

Instrument sets have been validated with Solus products only. **No additional items should be added to the set for sterilization.**

Prior to sterilization, loosen screwtop on Inserter to allow for proper steam penetration.

Method	Cycle Type	Minimum Temperature	Exposure Time	Minimum Drying Time
Steam	Pre-vacuum	270°F (132°C)	4 minutes	30 minutes

For surgical technique manual or additional information regarding the Alphatec Solus Anterior Lumbar Interbody Fusion (ALIF) Spinal Spacer System, please contact Alphatec Spine, Inc. Customer Service at +1 (800) 922-1356 or +1 (760) 431-9286.

Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician.

The Alphatec Solus™ product is not sponsored, endorsed by, or in any way affiliated with SOLAS® or The Society of Lateral Access Surgery.

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