

SAFEOP ORDERING DETAILS

Set Description	Set ID	Required/Optional
SafeOp Capital Equipment Set	SAFEOP	REQUIRED. CONTAINS AIX TABLET, SAFEOP HEAD UNIT, POWER SUPPLIES, AND DATA CABLE
SafeOp Disposables Soft Case	SODISP	REQUIRED. CONTAINS NEEDLE ELECTRODE PACKS, SAFEOP CLIPS, AND SAFEOP PROBES.

SafeOp Electrode Packs

■ AIX1201A	SafeOp EMG/NMJ Needle Electrode Procedure Kit
AIX1202A	SafeOp EMG/NMJ Surface Electrode Procedure Kit
■ AIX1203A	SafeOp EMG/NMJ Needle Electrode & SSEP Needle Electrode Procedure Kit
AIX1204A	SafeOp EMG/NMJ Needle Electrode & SSEP Surface Electrode Procedure Kit
AIX1205A	SafeOp EMG/NMJ Surface Electrode & SSEP Surface Electrode Procedure Kit

SafeOp Stimulating Accessories

■ AIX1310-S	SafeOp Stimulating Clip, Sterile
■ AIX1320-S	SafeOp Stimulating Ball-Tip Probe, Sterile
AIX1330-S	SafeOp Insulated Dilator Kit, Sterile

■ Standard in SODISP Set

DISPOSABLE USE EXAMPLES BY CASE TYPE

Case Type	Electrode Pack	Stimulating Accessories
Screw Testing	AIX1201A or AIX1202A	AIX1310-S and/or AIX1320-S
Lateral Case: EMG	AIX1201A or AIX1202A	AIX1310-S, AIX1320-S, and AIX1330-S
Lateral Case: EMG & SSEP	AIX1203A, or AIX1204A, or AIX1205A	AIX1310-S, AIX1320-S, and AIX1330-S